

BERPU[®] Syringe

2-Part Syringe

Defined position of the plunger at volume "0" to feel when plunger is completely inserted

Smooth flow, tight, high transparency barrel

Defined safety stop to avoid accidental pull-out of plunger




Well-defined markings

CE 0123 STERILE EO

• Latex free, non-pyrogenic




No rubber, styrene or DEHP, latex- and silicone-oil-free



2-Part Syringe

Size	1ml, 2ml, 3ml, 5ml, 10ml, 20ml	
Material	Barrel	medical grade polypropylene
	Plunger	medical grade polypropylene
	Lubricant	Oleic acid amide
	Graduation	indelible ink
Type	Centric luer slip, eccentric Luer slip, Hypodermic needles in all standard sizes can be attached to the syringes.	
Intended Use	The Sterile Hypodermic Syringes for Single Use are intended to be used for medical purpose to inject fluid into or withdraw fluid from body.	
Latex free	Yes	
DEHP/PVC/BPA free	Yes	
Sterilization method	EO	
Shelf life	5 years	
Standard	ISO 7886:2017	
Certification	CE, ISO13485	

Technical Data Sheet



Unit Packaging:

Syringe uses a paper based top and co-extruded film bottom.

Packaging information for 2-Part Syringe

ITEM#	ML/CC	BOX QTY	CASE QTY	SIZE
S20101001	1ml	100/box	20/case	450x310x435
S20201001	2ml	100/box	12/case	385x270x452
S20301001	3ml	100/box	12/case	385x270x452
S20501001	5ml	100/box	12/case	562x294x456
S21001001	10ml	100/box	12/case	560x355x435
S22001001	20ml	50/box	12/case	552x395x535

BERPU®

Tel: 86 577 85501388

E-mail: sales@berpu.com

Website: www.berpu.com

Address: NO.14 Xingji Road, Yongxing Street, Longwan District,
325000 Wenzhou Zhejiang Province, P.R. China