

BERPU®


Sterile Pen Needles For Single Use



- Ultra-fine needle
-special cut for gentle injection
- Special silicon coating
-for advanced injection comfort
- Thinnest 34G shortest 4MM available.
- Compatible with all major pens in the market.



RW, TW, and ETW available



Larger and more rapid insulin flow


STERILE EO

CE 0123

2

Technical Data Sheet

BERPU®



Pen Needles

Size	29G to 34G, needle length from 4mm to 13mm	
Material	Hub	Medical grade PP
	Needle protector	Medical grade PP
	Hub protector	Medical grade PP
	Cannula	SUS304
	Lubricant	Silicon oil
Intended use	The pen needles for single use are used with insulin pen for subcutaneous injection of insulin.	
Latex free	Yes	
DEHP/PVD/BPA free	Yes	
Sterilization method	EO	
Shelf life	5 years	
Standard	ISO 11608-2	
Certification	CE, ISO13485, FDA	

