

BERPU®

Blood Collection Set

- Wide, flexible wings for better control.
- Different color codes for easy identification.

- Diverse combinations available.



Sharp & smooth cannula edges make penetration painless




Compatible with all available holders.

CE 0123

STERILE EO

Technical Data Sheet



1. Rubber sleeve 2. Non-patient end needle 3. Needle hub 4. Tubing
5. Wing hub 6. Patient end needle 7. Needle protector

Blood Collection Set

| | | |
|----------------------|--|----------------------|
| Gauge | 18G-27G | |
| Material | Wing Hub | Medical grade PVC |
| | Needle Protector | Medical grade PP |
| | Tubing | Medical grade PVC |
| | Cannula | SUS304 |
| | Luer adaptor hub | Medical grade ABS |
| | Rubber sleeve | Medical grade rubber |
| Intended use | The blood collection set is applicable for blood collection from patients. | |
| Latex free | Customized | |
| Sterilization method | EO | |
| Shelf life | 5 years | |
| Standard | ISO7864:2016 | |
| Certification | CE, ISO13485 | |

