

BERPU®


Safety Blood Collection Set



- Wide, flexible wings for better control
- Protective shield prevents needlesticks
- Activate safety shield with one hand
- Transparent tubing to view blood flashback
- Compatible with all available holders

STERILE EO

CE 0123



- 1.Needle protector 2-Patient end needle 3.Wing hub
4.safety sheath 5.Tubing 6.Luer connector 7.Luer adapter hub
8.Luer adaptor cannula 9.Rubber sleeve 10.Luer adaptor cap

Safety Blood Collection Set

Gauge	18G-27G	
Material	Wing Hub	Medical grade PVC
	Needle Protector	Medical grade PP
	Tubing	Medical grade PVC
	Cannula	SUS304
	Luer adaptor hub	Medical grade ABS
	Rubber sleeve	Medical grade rubber
Intended use	The safety blood collection set for single use are applicable for blood collection from patients.	
Latex free	Customized	
Sterilization method	EO	
Shelf life	5 years	
Standard	ISO7864:2016	
Certification	CE, ISO13485	

