

BERPU[®]


Medical Rubber Sleeve For Blood Collection Needles

Precision design
avoids bias



Leakage proof for multiple
puncture

Anti-coring




High elastic recovery

STERILE | EO

CE 0123

Technical Data Sheet



Medical Rubber Sleeve For Blood Collection Needles

Material	Natural Rubber/ Latex Free Rubber
Intended use	Part of blood collection needle
DEHP/PVD/BPA free	Yes
Shelf life	5 years

Unit Packaging: PE bag

Packaging information for Rubber Sleeve

Volume/pcs	Carton packaging size
100,000	28.5cm*29cm*34cm

BERPU®

Tel: 86 577 85501388

E-mail: sales@berpu.com

Website: www.berpu.com

Address: NO.14 Xingji Road, Yongxing Street, Longwan District,
325000 Wenzhou Zhejiang Province, P.R. China