


## Cannula



Ultra sharp tips  
-for smooth penetration  
with minimal trauma

Fully automatic grinding  
production line  
-ensures high consistency and  
precision of every cannula

High quality austenitic stainless  
steel as raw material for production

RW, TW & ETW are all available



Anesthesia Cannula  
(Extradural Cannula)



Anesthesia Cannula



Indwelling Cannula




Biopsy Cannula




Biopsy Cannula



Pencil Point Cannula




Back Hole Cannula



Special Cannula



Double Hole Cannula



Bone Marrow Biopsy Cannula





Endoscopic Cannula



Different specifications  
can be customized

Special cannulas for  
customizations

# Technical Data Sheet



## Cannula

Gauge	11G - 34G, 15.5 - 44.5mm
Material	Medical grade stainless steel SUS304
Surface treatment	Electro-polished
Intended use	Part of all types of needle
Shelf life	5 years
Standard	ISO9626:2016
Certification	CE, ISO13485, FDA

### **Packaging:**

#### Packaging information for Cannulas

CASE MATERIAL	DIMENSION(mm)
Plywood case	350x245x270

# BERPU®

Tel: 86 577 85501388

E-mail: sales@berpu.com

Website: www.berpu.com

Address: NO.14 Xingji Road, Yongxing Street, Longwan District,  
325000 Wenzhou Zhejiang Province, P.R. China