


BERPU®

Fine Hypodermic Needle



- Smooth needle tip
-reduce tissue damage
- High-quality stainless steel
- Two types structure:
-type 1 hypdermic needle
-type 2 capsule

Color marked hub
-for easy identification of needle size

STERILE | EO

CE 0123



Fine Hypodermic Needles

Size	27G,30G, 31G, 32G, 33G, 34G	
Material	Needle Hub	Medical grade PP
	Needle protector	Medical grade PP
	Hub protector	Medical grade PP
	Cannula	SUS304
	Lubricant	Silicon oil
Intended use	It is applicable for injection and extraction of fluid in human skin, subcutaneous muscles, veins, etc.	
Latex free	Yes	
Sterilization method	EO	
Shelf life	5 years	
Standard	EN ISO 7864:2016	
Certification	CE, ISO13485	

