


Cannula



Ultra sharp tips
-for smooth penetration
with minimal trauma

Fully automatic grinding
production line
-ensures high consistency and
precision of every cannula

High quality austenitic stainless
steel as raw material for production

RW, TW & ETW are all available



Anesthesia Cannula
(Extradural Cannula)



Anesthesia Cannula



Indwelling Cannula



Biopsy Cannula



● Different specifications
can be customized




Biopsy Cannula



Pencil Point Cannula



Back Hole Cannula




Special Cannula



● Special cannulas for
customizations



Double Hole Cannula



Bone Marrow Biopsy Cannula



Endoscopic Cannula



Technical
Data
Sheet



RB

LB

SB

Cannula

Gauge	11G - 34G, 15.5 - 44.5mm
Material	Medical grade stainless steel SUS304
Surface treatment	Electro-polished
Intended use	Part of all types of needle
Shelf life	5 years
Standard	ISO9626:2016
Certification	CE, ISO13485, FDA