

## *Dental Needles For Single Use*


Smooth surface with light silicone coating  
-minimal pain upon puncture

Universal plastic hub  
-metric standard and British  
standard threads, fit most syringes

Lancet bevel point at  
cartridge end  
-prevents blockage of  
anesthetic

Triangle bevel mark  
-for optimal injection position

# Technical Data Sheet



## *Dental Needles For Single Use*

Size	30G, 27G, 25G	
Material	Needle tube	SUS304
	Needle hub	Medical grade PP
	Needle sheath	Medical grade PP
Intended use	The dental needles for single use is used with dental cartridge syringes for injection of dental local anaesthetics.	
Latex free	Yes	
DEHP/PVD/BPA free	Yes	
Sterilization method	EO	
Shelf life	5 years	
Standard	ISO 7885-2010	
Certification	CE, ISO13485, FDA	