

BERPU®

Vacuum Blood Collection Tube

Sterile, smooth, high transparency, complete seal



Additive uniformly on tube inner surface





Color-coded closures for accurate tube selection, easy identification

Hot runner mold production, ensure product higher stability

● Accessories like tubes and closures can be sold separately according to customer requests

CE 0123

STERILE EO



Blood Collection Tube

Color	Red/Orange/Yellow/Green/wathet blue/Black /Light grey/light Red/light green/Lavender
Material	Medical grade PET,Medical grade rubber, Medical grade PE
Draw Volume (ml)	1,2,3,4,5,6,7,8,9,10,11,12,15,16,8.5
Addictives	No additive / Clot Activator / Gel& Clot Activator / Heparin Lithium / Heparin Sodium / Sodium citrate1:9 / SodiumCitrate1:4 / Sodium fluoride / EDTA Na ₂ ·2H ₂ O / EDTA Na ₂
Intended use	Venous blood sampling for clinical use
Sterilization method	E-beam sterilization
Shelf life	1 year
Standard	EN14820:2004