

BERPU[®]

Sterile Safety Pen Needles

For Single Use



Simple and
Beautiful appearance

- Ultra-fine needle
-special cut for gentle injection
- Large barrel
-easy handing and attachment
- Compatible with all major pens in the market.



Special silicon coating for
advanced injection comfort.

Patented automatic safety mechanism
-prevent needlestick

Before use

After use



CE 0123



STERILE EO



Safety Pen Needles

Size	29G to 31G, needle length from 4mm to 12mm	
Material	Needle protector	PP
	Needle base	PP
	Barrel	PP
	Spring	SUS304
	Cannula	SUS304
	Lubricant	Silicon oil
Intended use	The Safety Pen Needles for single use are used with insulin pen for subcutaneous injection of insulin.	
Latex free	Yes	
DEHP/PVD/BPA free	Yes	
Sterilization method	EO	
Shelf life	5 years	
Standard	ISO 11608-2, EN ISO23908:2011	
Certification	CE, ISO13485, FDA	