

BERPU[®] Syringe

3-Part Syringe


1ml / Low Dead Space



Specially formulated gaskets
and silicone lubrication
– assures smooth and
accurate motion

High-transparency
syringe barrel

Defined safety stop to
avoid accidental
pull-out of plunger



Well-defined markings



CE 0123

STERILE EO

• Latex free, non-pyrogenic



3-Part Syringe

Size	1ml,2ml,3ml,5ml,10ml,20ml,30ml,50ml,60ml	
Material	Barrel	medical grade polypropylene
	Plunger	medical grade polypropylene
	Gasket	latex free nature rubber
	Lubricant	Silicon
	Graduation	indelible ink
Type	Centric luer slip, Eccentric Luer slip, Luer Lock, Hypodermic needles in all standard sizes can be attached to the syringes.	
Intended Use	The Sterile Hypodermic Syringes for Single Use are intended to be used for medical purpose to inject fluid into or withdraw fluid from body.	
Latex free	Yes	
DEHP/PVC/BPA free	Yes	
Sterilization method	EO	
Shelf life	5 years	
Standard	ISO 7886:2017	
Certification	CE, ISO13485, FDA 510K	