


# BERPU<sup>®</sup> Syringe

## 2-Part Syringe

Defined position of the plunger at volume "0" to feel when plunger is completely inserted

Smooth flow, tight, high transparency barrel

Defined safety stop to avoid accidental pull-out of plunger




Well-defined markings


CE 0123

STERILE EO

• Latex free, non-pyrogenic



No rubber, styrene or DEHP,  
latex- and silicone-oil-free



## 2-Part Syringe

Size	1ml, 2ml, 3ml, 5ml, 10ml, 20ml	
Material	Barrel	medical grade polypropylene
	Plunger	medical grade polypropylene
	Lubricant	Oleic acid amide
	Graduation	indelible ink
Type	Centric luer slip, eccentric Luer slip, Hypodermic needles in all standard sizes can be attached to the syringes.	
Intended Use	The Sterile Hypodermic Syringes for Single Use are intended to be used for medical purpose to inject fluid into or withdraw fluid from body.	
Latex free	Yes	
DEHP/PVC/BPA free	Yes	
Sterilization method	EO	
Shelf life	5 years	
Standard	ISO 7886:2017	
Certification	CE, ISO13485	