



BERPU®

Insulin Syringes For Single Use



CE 0123 STERILE EO

• Latex free rubber



Insulin Syringes

Size	0.3ml, 0.5ml, 1ml	
Material	Plunger	Medical grade PP
	Hub	Medical grade PP
	Hub protector	Medical grade PP
	Needle protector	Medical grade PP
	Rubber piston	Latex-free rubber
	Cannula	SUS304
	Lubricant	Silicon oil
Type/scale	U40 and U100, 0.3ml, 0.5ml, 1ml	
Intended use	The Insulin Syringe for single use is used for subcutaneous injection of insulin.	
Latex free	Yes	
DEHP/PVD/BPA free	Yes	
Sterilization method	EO	
Shelf life	5 years	
Standard	ISO 8537:2016	
Certification	CE, ISO13485, FDA	